

# KPL Alkaline Phosphatase Labeled Streptavidin (AP-SA)

*molecular biology grade*

Catalog No.  
5950-0005 (475-3000)

Size  
1 mL

## DESCRIPTION

Streptavidin is a 60,000 dalton protein isolated from the bacterium *Streptomyces avidinii*. The use of streptavidin, rather than egg white avidin, as the bridging reagent ensures that these products demonstrate sensitivity, high specificity and low background.

Streptavidin has been shown to bind four molecules of biotin with high affinity ( $K_p=10^{-15} M^{-1}$ ). Electrophoretically pure streptavidin is linked to alkaline phosphatase using a modified proprietary method adapted from Voller, et. al. <sup>(1)</sup>.

## FORM

The conjugate is provided in a liquid buffer containing a proprietary stabilizer and an anti-bacterial agent. It is prepared with molecular biology grade chemicals and dispensed into RNase/DNase-free sterile vials.

## STORAGE/STABILITY

Store at 2-8°C. Stable for a minimum of 1 year from date of receipt at 2-8°C as an undiluted liquid. Dilute only immediately before use.

## ENZYME:PROTEIN RATIO

Molar phosphatase:streptavidin ratio = 1.34:1.

## APPLICATIONS

Phosphatase labeled streptavidin is suitable for use in Northern blotting, Southern blotting, plaque and colony hybridizations, *in situ* hybridization, and immunohistochemistry. This conjugate may also be used for ELISA and immunoblotting procedures (See References 2-5).

## SUGGESTED WORKING DILUTIONS

Different assay conditions require that serial dilutions of all reagents be performed to determine optimal working concentrations. Prepare the working dilution immediately before use. Storage at a working dilution may result in enzyme inactivation and performance loss. Do not use PBS as a diluent. For suggested starting dilutions, see the appropriate protocol.

## SUGGESTED PROTOCOLS

All steps are at room temperature unless otherwise noted

### Northern Blotting, Southern Blotting, Plaque and Colony Hybridizations:

Following hybridization with a biotinylated probe and post-hybridization washing:

1. Place membrane in a small container or hybridization bag, block with KPL Detector Block, or other appropriate blocking solution, for 30 minutes.
2. Dilute KPL AP-SA in blocking solution (1:500-1:1000 for chromogenic detection and 1:10,000-1:50,000 for chemiluminescent detection).
3. Incubate membrane for 20 minutes with diluted AP-SA. Use at least 0.25 ml per cm<sup>2</sup> membrane. The optimal dilution of AP-SA must be determined experimentally.
4. Transfer membrane to a clean container and wash with KPL Biotin Wash Solution. Wash 3 times for 5 minutes each using at least 0.4 mL wash solution per cm<sup>2</sup> membrane.
5. Detect using KPL BCIP/NBT, CDP-Star<sup>®</sup> Chemiluminescent Substrate, or other appropriate phosphatase substrate following appropriate protocols.

### In Situ Hybridization:

Following hybridization of tissue or cells with a biotinylated probe and post-hybridization washing:

1. Dilute KPL AP-SA 1:20-1:200 in an appropriate diluent and apply approximately 100 µL to specimen. Cover slide to prevent evaporation. Incubate in a 37°C humidified chamber for 20 minutes. The optimal dilution of KPL AP-SA must be determined experimentally.
2. Immerse slide in a coplin jar containing KPL Biotin Wash, or other appropriate wash solution. Wash 3 times for 5 minutes each.
3. Detect using KPL HistoMark<sup>®</sup> RED, KPL HistoMark<sup>®</sup> BLUE or other appropriate phosphatase substrate following manufacturer's instructions.
4. Counterstain, if desired, and mount slide.

# KPL Alkaline Phosphatase Labeled Streptavidin (AP-SA)

*molecular biology grade*

Catalog No.  
5950-0005 (475-3000)

Size  
1 mL

## Immunohistochemistry:

Following incubation of the specimen with primary and biotin-labeled secondary antibody:

1. Dilute KPL AP-SA 1:10-1:100 in KPL Goat Serum Block or 10% normal goat serum in 0.01 M Tris, pH 7.65. Flood the slide with diluted KPL AP-SA.
2. Incubate at room temperature for 30 minutes. The optimal dilution of KPL AP-SA must be determined experimentally.
3. Soak the slide in Tris-HCl for 5 minutes.
4. Detect using KPL HistoMark® RED, KPL HistoMark®BLUE or other appropriate phosphatase substrate.
5. Counterstain if desired and mount.

## PRODUCT SAFETY AND HANDLING

This product is considered non-hazardous as defined by The Hazard Communication Standard (29 CFR 1910.1200). Avoid contact with skin and eyes. In case of contact or spillage, clean with copious amounts of water. Product may be disposed via a sanitary sewer.

## REFERENCES

1. Voller, A., et. al. (1976). A Manual of Clinical Immunology. American Society for Microbiology.
2. Brigati, D.J., et. al. (1983). *Virology*, 126: 32-50.
3. Harlowe, E. and Lane, D. (1988). Antibodies: A Laboratory Manual. Cold Spring Harbor Laboratory Press, NY.
4. Sambrook, J. et. al. (1989). Molecular Cloning: A Laboratory Manual. Cold Spring Harbor Laboratory Press, NY
5. Crowther, J.R. (1995). Methods in Molecular Biology: ELISA: Theory and Practice. Vol. 42, Humana Press, NJ.

## RELATED PRODUCTS

## CAT. NO.

KPL Detector Block	5920-0004 (71-83-00)
KPL Biotin Wash	5960-0015 (50-63-06)
KPL Goat Serum Block	5560-0007 (71-00-27)
KPL BCIP/NBT Substrate	5420-0038 (50-81-18)
KPL HistoMark Red	5510-0036 (55-69-00)
KPL HistoMark Blue	5510-0037 (55-70-00)

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis.